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Instruction Manual

Thank you for choosing our product. We hope you enjoy using the appliance.

Contents

General Notes	35
Symbols in these Instructions for Use	35
Children.....	35
Special Safety Instructions For This Machine	35
Handling Batteries	36
Inserting the Batteries	36
Location of Controls	36
Technical Specifications	36
Information about the Analysis Scale	37
Use	37
Weighing Function	37
Saving Personal Data	37
Measurement and Analysis.....	37
Assessment of the Results	38
Body Fat.....	38
The Importance of Water for our Health.....	38
Guideline values for body fat and body water content ..	38
BMR (Basal Metabolic Rate / Energy expenditure).....	39
Reference values of the daily intake of	
kilocalories without physical activities	39
BMI (Body Mass Index)	39
Reference values for the BMI.....	39
Display Error Messages	39
Cleaning	39
Disposal	39
Disposal of Batteries.....	39
Meaning of the "Dustbin" Symbol	40

General Notes

Read the operating instructions carefully before putting the appliance into operation and keep the instructions including the warranty, the receipt and, if possible, the box with the internal packing. If you give this device to other people, please also pass on the operating instructions.

- The device is designed exclusively for private and not for commercial use.
- Only use the device for the intended purpose.
- Do not use a damaged device.
- Keep away from heat, direct sunlight, damp (never immerse in liquids) and sharp edges.
- Do not try to repair the appliance on your own. Always contact an authorized technician.
- Pay careful attention to the following "Special Safety Instructions".

Symbols in these Instructions for Use

Important information for your safety is specially marked. It is essential to comply with these instructions in order to avoid accidents and prevent damage to the machine:

WARNING:

This warns you of dangers to your health and indicates possible injury risks.

CAUTION:

This refers to possible hazards to the machine or other objects.

NOTE:

This highlights tips and information.

Children

- In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach.

WARNING!

Do not allow small children to play with the foil. There is a **danger of suffocation!**

- Children should be supervised at all times in order to ensure that they do not play with the device.

Special Safety Instructions For This Machine

WARNING:

If the glass surface is moist or covered with water droplets, please wipe the surface dry as there is a **danger of slipping**.

CAUTION:

- Handle the balance with care, as it is a precise measuring instrument. Do not drop it or jump on it.
- Do not load the appliance with more than 180 kg (396.8 lb, 28: 3½ st:lb) to avoid damaging the device.
- If the scales are carried or not going to be used, please store them only in a horizontal position.
- The body fat analysis scale is not suitable for the following persons:
 - Children under 6 years and adults over 100 years
 - Persons with symptoms of fever, oedema, or osteoporosis
 - Persons undergoing treatment for dialysis
 - Persons who use a pacemaker
 - Persons who take cardiovascular medication
 - Pregnant women
 - Sportspersons who perform over 5 hours of intensive sporting activities each day
 - Persons with a heartbeat lower than 60 beats/minute

36 English

- Persons who are taking vaso-dilating or vaso-constrictive medications
- Persons with considerable anatomic deviations of the legs in relation to total body size (length of the legs considerably shortened or lengthened).

Handling batteries

⚠ WARNING:

- Do not expose the batteries to high temperatures or direct sunlight.
- Never throw batteries into the fire. There is a **danger of explosion!**
- Keep batteries away from children. They are **not** a toy.
- Batteries may be swallowed, which can be life-threatening. If a battery has been swallowed, seek medical help immediately.
- Do not open the batteries by force.
- Avoid contact with metallic objects. (Rings, nails, screws et cetera.) There is a **danger of short-circuiting!**
- As a result of a short-circuit batteries may heat up considerably or even catch fire. This may result in burns.
- For your own safety the battery poles should be covered with adhesive strips during transport.
- If a battery leaks do not rub the liquid into your eyes or mucous membranes. If hands come into contact with the liquid they should be washed thoroughly and eyes rinsed out with clear water; if symptoms persist please consult a doctor.

Inserting the Batteries

📌 NOTE:

Where required, the supplied batteries have been protected for transport with plastic foil.

1. Remove this foil before inserting the batteries into the appliance.
2. Open the battery compartment on the back.
3. Insert 2 batteries of type AAA/R03. Please ensure that the polarity is correct (see the floor of the battery compartment!)
4. Close the battery compartment again.
5. If the device is not used for a prolonged period of time, please take the batteries out in order to prevent leakage of battery acid.
6. If the message "Lo" appears on the display, replace the batteries.

⚠ CAUTION:

- Different battery types or new and used batteries must not be used together.
- Batteries are not to be disposed of together with domestic wastes. Dispose of the used batteries at the relevant collection bin or at your dealer's.

Location of Controls

- 1 UP-button ▲
- 2 SET-button ■
- 3 DOWN-button ▼
- 4 Display
- 5 Contact Surfaces

Display Indications	Meaning
kg (lb / st:lb)	Body weight
P0 - P9	User memory space
	Male
	Female
cm (ft:in)	Body height
age	Age
fat	Body fat display
TBW	Body water display
	Muscle mass display
	Bone weight display
%	Percentage – body fat / body water / muscle mass
Underfat Healthy Overfat Obese	Rating: underweight – normal weight – overweight - obese
KCAL	Basal metabolic rate of the daily energy requirement in kilocalories (BMR)
BMI	Body Mass Index

Technical Specifications

Model:PW 5644 FA
 Power supply: 2 Batteries 1.5 V Type "AAA" / "R03"
 Net weight of device: approx. 1.57 kg
 Memory capacity: 10
 Mass alignment: from 80 to 220 cm (2'7½" – 7'2½" ft:in)
 Age: from 6 to 100 years
 Minimum weight: from 5.0 kg (11.0 lb / 0: 11 st:lb)
 Maximum weight: 180 kg (396.8 lb / 28: 3½ st:lb)
 Weight measurement unit: .. 100 g (0.2 lb / ¼ st:lb) increments
 Body fat measurement unit: 0.1 %
 Body fluid graduation: 0.1 %
 Recommended operating environment:
 Temperature: 0 °C – 40 °C / 32 °F – 104 °F
 Relative air humidity: < 85 %

The right to make technical and design modifications in the course of continuous product development remains reserved.

This device has been tested according to all relevant current CE guidelines, such as electromagnetic compatibility and low voltage directives, and has been constructed in accordance with the latest safety regulations.

Information about the Analysis Scale

WARNING:

Do not use the scales if you use a pacemaker or another medical implant with electronic components. The function of the implant may be impaired by the measurement.

The main principle of this device is based on the measurement of the electric impedance within the human body. In order to determine the actual body data, the device leads a weak, unnoticeable electric signal through your body.

The signal is generally safe! However, please read our special safety instructions!

The method of measurement is referred to as "bioelectric impedance analysis" (BIA), a factor which is based on the relationship between body weight and associated body water percentages as well as on other biological data (age, gender, body weight).

Use

NOTE:

- You will find the buttons for the operation on the top edge of the scale.
- The scales are factory set with the units "kg" and "cm".

To change the settings to "lb" or "st:lb" as well as "ft:in", proceed as follows:

- Switch the scales on by applying gentle pressure on the glass surface with your foot or toe.
- Press the UP button ▲, until the required unit of weight is displayed. The unit for the height also automatically changes.

Weighing Function

- Place the scale on a hard, flat surface (avoid carpets). An uneven floor contributes to inaccurate measurements.
- Stand on the scales with both feet. The device switches on automatically.
- Distribute your weight equally and wait for the scale to determine your weight.
- The digits on the display will flash before your exact weight will be shown.
- In order to extend the battery life, the device will automatically turn off 15 seconds after you have stepped off the scale.
- If there is no weight determination while "0.0" is shown, the scale will shut off automatically after 15 seconds.

Saving Personal Data

NOTE:

- The scale has default values that you can change, as described below.
- Carry out the following steps in sequence and within a few seconds. Otherwise, the device will automatically cancel input mode.
- Press the button SET. The memory space (P0) will flash on the display. Confirm the setting by pressing SET or select another memory space by pressing the UP / DOWN button. Then press SET to go to the next setting.
- The symbol for "male" or "female" will flash. Confirm the setting by pressing SET or use the UP / DOWN button to choose the symbol for your gender. Then press SET to go to the next setting.
- The height setting will flash. Select your height by pressing UP or DOWN button. Holding down one of the buttons will accelerate the process. Confirm your entry by pressing SET.
- The setting AGE will flash. Select your age by pressing UP or DOWN button. Holding down one of the buttons will accelerate the process. Confirm your entry by pressing SET.
- The short display " - - - " will terminate the entry mode for this memory location. Then you will see the memory location, the symbol of the gender and the weight "0.0 kg" (0.0 lb / 0: 0 st:lb) in the display.

You now have the following options:

- You can change the units for weight and type. Press the UP button ▲, until the desired unit of weight is displayed. The unit for the height automatically changes.
- Press the SET button to see the selected settings again.
- Allocate further memory spaces. There are 10 spaces available in total. Press SET and change the space by pressing the UP or DOWN button. Proceed further as described above.
- If you want to make changes, press SET (the memory location will blink). Select the stored position using the UP or DOWN button. Proceed further as described above.
- You can immediately start measuring.

Measurement and Analysis

Record your personal data in the memory.

NOTE:

- The analysis scale will only work correctly when you are standing on the contact surfaces barefoot and with slightly moistened foot soles, if possible. Completely dry foot soles can lead to unsatisfactory results as they have a conductivity that is too low.
- If necessary wipe the contact surfaces to improve results.

i NOTE:

- Due to the significant fluctuation of body fat, it is difficult to attain a standard value. On average the optimal testing time is between 6 p.m. and 8 p.m., since the body fat remains very consistent during this period of time.
- Press the SET button. Use the UP / DOWN button to select the stored position at which your data is saved.
- Your preset data briefly appears in the display. Wait until the display comes to a standstill: memory location, the symbol of the gender and the weight "0.0 kg" (0.0 lb / 0: 0 st:lb).
- Stand on the contact surfaces (5) of the scales with both feet.
- Firstly your weight will flash on the display.
- You'll see a treadmill "o o o o" in the display during the analysis. After a few seconds, the following results will be shown:
 1. The body fat proportion in % (**fat**) and the rating
 2. The proportion of body fluid % (**TBW**)
 3. The proportion of muscle mass in % (symbol )
 4. The bone weight in kg (lb / st:lb) (symbol )
 5. The basal metabolic rate of the daily energy requirement in **KCAL**
 6. The Body Mass Index (**BMI**)
- Afterwards the determined data will be shown repeatedly. Then the scale will automatically turn off.

Assessment of the Results

There are heavy persons with a well-balanced percentage between fat and muscle mass, and there are seemingly slim persons with significantly high body fat values. Fat is a factor that puts serious strain on the body, in particular the cardiovascular system. Therefore, controlling both values (fat and weight) is a prerequisite for health and fitness.

Body Fat

Body fat is an important component of the body. From protecting the joints as a vitamin repository to regulating the body temperature, fat plays a vital role in the body. Thus the goal is not to radically reduce the proportion of body fat, but rather to achieve and maintain a healthy balance between muscle and fat.

The Importance of Water for our Health

Approximately fifty-five to sixty percent of the human body consist of water, depending on age and gender in varying amounts.

Water has to fulfill several tasks within the body:

- It is the building block of our cells. That means that all the cells in our body, whether skin cells, gland cells, muscle cells, brain cells or others, are functional only if they contain enough water.
- It is a solvent for important substances in our body.
- It is a means of conveyance for nutrients, the body's own substances and metabolic products.

Reference values for your fitness profile can be seen in the benchmark table. It shows the dependency of body fat data on age and gender. Always pay attention to your body's water balance, which, as well as the fat level, undergoes strong fluctuations throughout the day. Therefore, to obtain consistent data, always measure at the same time of day.

The calculated results of the body fat and water level should only be considered as reference values. They are not recommended to be used for medical purposes. Please always consult a doctor if you have questions about a medical treatment or a diet in order to achieve a normal body weight.

Guideline values for body fat and body water content as a percentage (approximate figures)

Age (Years)	Women		Men		
	Body fat	Body water	Body fat	Body water	
< = 30	4.0 – 20.5	66.0 – 54.7	4.0 – 15.5	66.0 – 58.1	Underweight
	20.6 – 25.0	54.6 – 51.6	15.6 – 20.0	58.0 – 55.0	Normal
	25.1 – 30.5	51.5 – 47.8	20.1 – 24.5	54.9 – 51.9	Overweight
	30.6 – 50.0	47.7 – 34.4	24.6 – 50.0	51.8 – 34.4	Obese
> 30	4.0 – 25.0	66.0 – 51.6	4.0 – 19.5	66.0 – 55.4	Underweight
	25.1 – 30.0	51.5 – 48.1	19.6 – 24.0	55.3 – 52.3	Normal
	30.1 – 35.0	48.0 – 44.7	24.1 – 28.5	52.2 – 49.2	Overweight
	35.1 – 50.0	44.6 – 34.4	28.6 – 50.0	49.1 – 34.4	Obese

BMR (Basal Metabolic Rate / Energy expenditure)

The Basal Metabolic Rate (BMR) characterizes the energy expenditure. The energy expenditure is used to describe the amount of energy in in KCAL that the body requires per day during complete rest, at an indifference temperature (28-30 °C) and on an empty stomach for the maintenance of its function (e.g. while sleeping).

Reference values of the daily intake of kilocalories without physical activities

Age (Years)	Women		Men	
	Weight (kg)	KCAL	Weight (kg)	KCAL
6-17	50	1265	60	1620
18-29	55	1298	65	1560
30-49	60	1302	70	1561
50-69	60	1242	70	1505

BMI (Body Mass Index)

The Body Mass Index (BMI) refers the body weight to the square of the body height.

BMI = Weight in kg: (body height in m x body height in m)

Example:

A woman weighs 55 kg and is 1.68 m tall. Her BMI is calculated as follows:

$$55 \text{ kg} : (1.68 \text{ m} \times 1.68 \text{ m}) \approx 19.5 \text{ kg/m}^2 \rightarrow \text{Normal weight}$$

The value "square of the body height" is in no relation to the body surface. The BMI is merely a rough reference value, as it neither considers the build and gender, nor the individual composition of the body mass from the fat and muscle tissue of a human.

Reference values for the BMI

BMI	Under-weight	Normal	Over-weight	Obese
	< 18.5	18.5 - 25	25 - 30	> 30

Please note that the interpretation of the BMI indicates excess weight with very muscular bodies (bodybuilder). The reason for this is that the far above-average muscle mass is not considered in the BMI formula.

Display Error Messages

NOTE:

If an error message appears on the display, wait until the device is turned off. Turn it on again to restart the operation. If the error display has not disappeared, take the batteries out of the battery compartment for 10 seconds.

Display Indications	Meaning
	Self test (after inserting the batteries) for perfect functioning of the display
Lo	Battery is low
Err	The weight lies above the allowed maximum weight.
Err2	The percentage of body fat content lies above or below the limits of the measuring scale.

Cleaning

CAUTION:

- Do not immerse the balance in water: the electronic components can be damaged.
- Do not use a wire brush or any abrasive items.
- Do not use any acidic or abrasive detergents.
- Clean the device with a slightly moist cloth. Make sure that no humidity or waters seeps into the device!
- Use a commercially sold detergent or soap for cleaning.
- Dry the device well after cleaning with a dry and soft cloth.

Disposal

Disposal of Batteries

Your appliance contains batteries. If these batteries are spent, follow the instructions below:

CAUTION:

Batteries/rechargeable batteries must not be disposed of with household waste!

- The batteries must be removed from the appliance.
- Take spent batteries to the appropriate collection point or to a dealer.
- Your town or local authority can provide information about public collection points.

This symbol can be found on batteries/rechargeable batteries which contain hazardous substances:

- Pb = contains lead
- Cd = contains cadmium
- Hg = contains mercury
- Li = contains lithium

40 English

Meaning of the "Dustbin" Symbol

Protect our environment: do not dispose of electrical equipment in the domestic waste.

Please return any electrical equipment that you will no longer use to the collection points provided for their disposal.

This helps avoid the potential effects of incorrect disposal on the environment and human health.

This will contribute to the recycling and other forms of reutilisation of electrical and electronic equipment.

Information concerning where the equipment can be disposed of can be obtained from your local authority.

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